



**Complete line for clinical,
Rehabilitative use and sports medicine use**

RAM SRL
VIA GERMANIA 14 - 35127 PADOVA - ITALY
Tel. +39 049 870 33 55
bisellom@ram-treadmills.com

THE BEST CHOICE FOR YOUR PERFORMANCE

PERFECT FOR CLINICAL USE, REHABILITATION, AND SPORTS MEDICINE

RAM treadmills are designed to meet the rigorous demands of healthcare and sports environments. They offer a range of features aimed at ensuring optimal performance and maximum patient safety.

SUPERIOR COMFORT:

The triple-layer running belt is engineered to provide maximum energy absorption, making every workout as comfortable as possible.

ADVANCED SAFETY:

The non-slip rubber surface and parallel bars ensure excellent stability and grip, minimizing the risk of accidents.

UNMATCHED DURABILITY:

The reinforced frame guarantees long-lasting durability and reliable performance, even in high-usage environments.



CHOOSE YOUR MODEL

	870A	870S	870C
Leinght	205 cm	230 cm	252 cm
Width	80 cm	89 cm	96 cm
Height	112 cm	115 cm	115 cm
Running surface	50x150 cm	58x170 cm	65x190 cm
Weight	165 Kg	220 Kg	245 kg
Patient capacity	200 Kg	200 Kg	200 Kg
Speed range	0-20 Km/h	0-25 Km/h	0-32 km/h
Speed increments	0,1 Km/h	0,1 Km/h	0,1 Km/h
Incline range	0-25%	0-25%	0-22%
Incline increments	0,5%	0,5%	0,5%
Speed motor	Asynchronous three phase 2,0 HP, AC (nominal power)	Asynchronous three phase 3,0 HP, AC (nominal power)	Asynchronous three phase 4,0 HP, AC (nominal power)
Elevation motor	90Watt DC (nominal power)	90Watt DC (nominal power)	90Watt DC (nominal power)
Running belt	Impact absorbing (thickness 3,1mm)	Impact absorbing (thickness 3,1mm)	Impact absorbing (thickness 3,1mm)
Surface color	Embossed white	Embossed white	Embossed white
Handrails color	Black	Black	Black
Belt color	Black	Black	Black
Safety devices	2 red stop switches	2 red stop switches	2 red stop switches
Serial Port	RS232 asynchrone	RS232 asynchrone	RS232 asynchrone
Connector	DB9 F	DB9 F	DB9 F
Baude rate	4800 baud rate	4800 baud rate	4800 baud rate
Operation	Remote device interface	Remote device interface	Remote device interface
Packaging	Carton box on wood pallet	Carton box on wood pallet	Carton box on wood pallet

OPTIONAL



Short handrails



Display 15,6" touchscreen



Pediatric handrails
(disponible for 870A/870S)



Safety lanyard
Autostop function